

[Request a Quote](#)[◀ Previous page](#)

STS-UV Spectrometers


Please contact us for further details.

[Get in Touch](#)

Details

The STS-UV spectrometer (190-650 nm) is a versatile choice for applications ranging from low-concentration absorbance measurements to high intensity laser characterization. Its rugged design, compact size and great unit-to-unit reproducibility make STS attractive for integration into other devices and setups where a small footprint is desired.


- Cost-effective – optical resolution of 1.5 nm is comparable to more expensive setups
- Reproducible – attractive design for integration and scalability
- Multiple slit options -- STS-UV models are available with 10 μm , 25 μm , 50 μm , 100 μm or 200 μm slits



Using the Right UV Light Saves Lives >

In the battle against COVID-19, UV sterilization techniques can extend the life of high-demand personal protective equipment (PPE). But does your UV lamp provide enough power to disinfect face masks safely?

Spectrometer Model: Flame-S | Slit Size: 25 μm | Grating Type: 600



Specifications

Select model to see full specifications

Spectroscopy:

Wavelength Range:

190nm - 650nm

Optical Resolution:

1.0 nm
1.5 nm
12.0 nm
3.0 nm
6.0 nm

Integration Time:

10 μ s - 10s

Dark Noise:

\leq 3 counts RMS

Dynamic Range:

5 x 10⁹ (system
10 s max integration)
~4600 single acquisition

Input Fiber Connector:

SMA 905

Signal to Noise Ratio:

>1500:1 (full signal)

Stray light:

≤0.25%

Detector:**Detector:**

CMOS

Entrance slit:

- 10 µm
 - 100 µm
 - 200 µm
 - 25 µm
 - 50 µm
-

Detector Collection Lens:

No

Physical:**Dimensions:**

40 mm x 42 mm x 24 mm

Weight:

- ~ 60 g
- ~60 g

Show less



Related Products